



TheraTogs™ are a revolutionary orthotic undergarment and strapping system that gives clients with sensorimotor impairment an exciting new modality for improving postural alignment and stability, movement skill and precision, and joint stability.

**TheraTogs™ are designed and developed
by Beverly Cusick, PT, MS.**

TheraTogs™ provide a foam-lined, Latex®-free, breathable, comfortable, Velcro®-sensitive “external skin” for clinicians engaged in sensorimotor training of clients with developmental disabilities, and for clients in rehabilitation following brain injury, stroke, orthopedic surgery, or neurosurgery.



Use TheraTogs™ to:

- ✓ Promote optimum bone and joint development
- ✓ Provide joint stability and increased body awareness
- ✓ Improve posture, balance, gait, and movement skills
- ✓ Reduce the postural correction workload that is often expected of foot and ankle orthoses
- ✓ Improve motor performance during therapy sessions
- ✓ “Send the therapist’s hands home” in support of the therapy program between sessions

Custom fit each client with a cozy, non-slip TheraTogs™ set and an elastic strapping and support system designed to address that client’s specific musculoskeletal, postural, balance, or movement issues. Select from a variety of strap designs, lengths, widths, calibers of stiffness, and placement options to optimize their therapeutic and orthotic influences.

TheraTogs™ Indications:

- Ataxia
- Hypotonia
- Athetosis
- Weakness
- In-toed gait
- Out-toed gait
- Joint deviations
- Joint laxity
- Muscle recruitment problems
- Sensory integration issues
- Postural mal-alignment
- Muscle imbalance

TheraTogs™ Unique Features:

- Vertical resistance in the fabric promotes postural stability
- Horizontal fabric flexibility offers easy fit and client comfort – in all climates
- Cuffs extend strapping applications to lower and upper extremities
- Client enjoys toilet access without doffing TheraTogs
- Versatility by design:
- Custom sizing options are built into trim lines
- Hipster can be made unilateral
- Entirely adaptable strapping system
- Marker dots simplify caregiver’s re-application at home

You can choose from:

- ❑ **TheraTogs Full Body Kits** for patients with sensory issues, motor or balance impairment – any issue that requires strapping anchors to the full torso.
- ❑ **Clinical Sizing Kits** for clinics and therapy departments who want a permanent resource for designing TheraTogs solutions and testing size and effect for a range of patients.

You'll find detailed sections in this catalogue about these exciting items along with other information about TheraTogs applications.

TECHNICAL SPECIFICATIONS

TheraTogs™ can be worn directly on the skin, out of sight under the client's regular attire. When snugly applied, they comfortably adhere to the skin surface, which allows proper strapping applications to effectively influence the client's movement, stability, posture, and gait.

TheraTogs™ are made of a patented composite fabric consisting of nylon and spandex, with a foam backing made of an aqueous-based elastomeric urethane. TheraTogs are hand-washable. TheraTogs™ are Class I medical devices intended to be issued by, and applied under the supervision of, a licensed healthcare practitioner.

WARNING. THERATOGS ARE NOT MADE OF FLAME-RETARDANT MATERIALS, AND ARE NOT INTENDED TO BE WORN WHILE SLEEPING.



THERATOGS™ IS A TRADEMARK OF PROGRESSIVE GAITWAYS, LLC.
ALL RIGHTS RESERVED.
U.S. AND FOREIGN PATENTS PENDING.

For sales and further information, contact:

Rehab Logistics

Telephone: 00353 46 9027935 Fax: 00 353 46 9072065

Caroline Mooney: 00 353 87 9163457

caroline.mooney@rehabenterprises.ie