

stowaway

Positive Momentum™



Remploy are at the forefront of wheelchair design and development and as such are the choice for mobility. Remploy Mobility have a long established reputation for innovation in wheelchair design and construction.

We pride ourselves on our customer focus and after sales service. Our continuous programme of product development and marketing focus enables us to deliver innovative solutions to suit customer need; from the active user through to highly specialised clinical need. We are committed to the continued delivery of our **Positive Momentum™**.

ultra compact lightweight folding wheelchair

The Stowaway folding wheelchair allows the occupant to travel in style. Whether going on holiday, out for the day shopping or simply visiting friends this highly adaptable, ultra compact, lightweight (only 14kgs) wheelchair is ideal for the occasional user.

- ✦ Available in occupant and attendant propelled versions
- ✦ Unique 2 way folding system
- ✦ Folds down to the size of a large holdall
- ✦ Choice of 3 colours available
- ✦ 17" Seat Width
- ✦ New 'Fusion' padded upholstery for added comfort
- ✦ Fits 330 x 840 x 510 storage space, weighs less than 14kg
- ✦ Low maintenance steel frame
- ✦ Occupant weight limit 136kgs
- ✦ ISO9002 quality assured



Positive
Momentum



stowaway

technical information



Features

- ◆ Folds down to the size of a large holdall.
- ◆ Stowaway armrests provide side transfer by folding down to seat level.
- ◆ Available in Azzurri Blue, Midnight Blue and Silver.
- ◆ Stowaway footrests provide user access by folding up and swivelling sideways.
- ◆ Footrest from 510 - 660mm (Height Adjustable).



Weight (Kg)

Model Ref.	Total Weight	Weight for Lifting
AP	13.5	12.5
OP	14.5	10.0

Corridor Width Required to Turn (mm)

Model Ref.	Through 90° From a Door	Through 360° Between Walls
AP	990 (0.99m)	1580 (1.58m)
OP	990 (0.99m)	1180 (1.18m)

AP includes turning space for carer

Maximum Occupant Weight Tested to 136Kg (220lbs)

Model Ref: AP (Attendant Propelled) - OP (Occupant Propelled)

Seating Area Dimensions

Model Ref.	Seat Ref. (Inch)		Effective Measurements (mm)		
	Width	Depth	Width	Depth	Height
AP OP	17	16	465	440	415

Overall Dimensions – Chair Open and Folded

Model Ref.	Width Ref (inches)	Chair Open(mm)			Chair Folded(mm)		
		Width	Length	Height	Width	Length	Height
AP OP	17	595 595	915 940	940 940	330 330	840 840	510 510

Rear Wheel Information (mm)

Model Ref.	Fitting	Diameter	Width	Tyre	Handrim
AP	Fixed	315	40	Puncture Free	None
OP	Quick Release	510 (20")	31 25	Pneumatic Puncture Free	19mm Diameter MCP (carefree)

Remploy Healthcare reserve the right to change product specification without prior notice.



GANDON ENTERPRISES LTD
A PEARSON GROUP COMPANY

Phone: 0469027935

Fax: 0469072065

Remploy Healthcare

quality



service



innovation



www



sales

